



## Product Information

---

### **Rat Microsomal [triglyceride] Transfer Protein (MTP)**

Catalog No. RB-RMTP

**Description** Purified MTP from rat liver. Useful as a source of MTP for an activity assay.

**Source** Rat liver

**Unit Size** 0.5 ml at 1.7 mg/ml total protein

**Storage and Handling** Store at -80° C. Avoid freeze/thaw cycles.

**Caution:**

Wear appropriate gloves, protective clothing and eyewear.  
Follow safe laboratory practices.

### **Related Products**

RB-MTP Roar MTP Activity Assay Kit, 100 assays

RB-HMTP Roar Human Microsomal [triglyceride] Transfer Protein (MTP)

**For Research Use Only. Not for Diagnostic or Therapeutic Purposes.**

---

Roar Biomedical, Inc., Audubon Biomedical Center, 3960 Broadway, New York, NY 10032 USA  
Tel: +1 (212) 280-2983 • Fax: +1 (212) 280-2968 • info@roarbiomedical.com • www.roarbiomedical.com

©2016 Roar Biomedical, Inc. All rights reserved. This information is subject to change without notice.