




Product Information

Microsomal [triglyceride] Transfer Protein (MTP), partially purified

Catalog No. RB-HLV

Description	Partially purified MTP from human liver. Useful as a source of MTP for an activity assay.
Source	Human liver
Unit Size	0.5 ml at 20 mg/ml total protein
Storage and Handling	Store at -80° C. Avoid freeze/thaw cycles. Caution: Human Tissue Product HBsAg & HIV-1/2 antibody non-reactive. Treat all human tissue products as hazardous. Shelf life: 3 years, if stored properly. Wear appropriate gloves, protective clothing and eyewear. Follow safe laboratory practices.

Figure 1. MTP transfer activity using the Roar MTP Activity Assay Kit (cat. #RB-MTP) and partially purified human MTP (cat. #RB-HLV)



Related Product

RB-MTP Roar MTP Activity Assay Kit, 250 assays

For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

Roar Biomedical, Inc., Audubon Biomedical Center, 3960 Broadway, New York, NY 10032 USA
Tel: +1 (212) 280-2983 ▪ Fax: +1 (212) 280-2968 ▪ info@roarbiomedical.com ▪ www.roarbiomedical.com

©2018 Roar Biomedical, Inc. All rights reserved. This information is subject to change without notice.